




All-In-One Snare Probe

5 functions in one device

- Simple and user-friendly, seamless procedure integration
- Enables a wide range of procedures without instrument change
- Enables hybrid techniques, such as Hybrid-ESD, without altering the usual workflow



All-In-One Snare Probe

Details and components

Thanks to its innovative design, the All-in-One Snare Probe combines five essential functions for endoscopic procedures in the gastrointestinal tract in a single instrument:

- marking & coagulation
- incision
- dissection
- submucosal re-injection
- polypectomy

Precise control and the synchronously rotating snare ensure comfortable handling and procedure execution. By integrating a snare and a knife into a single instrument, this probe is optimally designed to meet the demands of modern endoscopic procedures.



Technical specifications

ALL-IN-ONE SNARE PROBE	Handle type	Snare shape	Snare opening width [mm]	Snare type	Working length [mm]	Tip length [mm]	Required working channel [mm]	Ref.No.
All-In-One Snare Probe	Rotational	Oval	25	Soft	2400	1.5	3.2	L24RA2425SF

Application



Marking/coagulation



Incision and dissection



Submucosal re-injection



Polypectomy

Source: Leo Medical

Ovesco Endoscopy AG
Friedrich-Miescher-Str. 9 | D-72076 Tuebingen
+ 49 (0) 7071 96528-160 | service@ovesco.com
www.ovesco.com

LeoMed
乐奥医疗
MAKE LIFE MORE HEALTHY

Leo Medical Co.,Ltd.
www.leomed.cn/en

